
MATERIAL SAFETY DATA SHEET
ACTICOAT PRODUCT RANGE
MSDS No. 179

SECTION 1 – IDENTIFICATION OF THE MATERIAL AND SUPPLIER

Product
Product Name ACTICOAT PRODUCT RANGE

Other Names

 ACTICOAT Antimicrobial Barrier Dressing
 ACTICOAT 7 Antimicrobial Barrier Dressing
 ACTICOAT ABSORBENT Antimicrobial Dressing
 ACTICOAT MOISTURE CONTROL Absorbent Antimicrobial Dressing
 ACTICOAT SITE Antimicrobial Barrier Dressing

Manufacturer's Product Code Use

 Various
 Antimicrobial wound dressing.

Supplier
Company Smith & Nephew Pty. Limited
Address 315 Ferntree Gully Road, Mount Waverley, Victoria 3149
Telephone No. Customer Service:
 (03) 8540 6777 (Toll Free 13 13 60)
Emergency Tel. Nos. (03) 8540 6777 (business hours)

SECTION 2 – HAZARDS IDENTIFICATION

Hazard Classification NON-HAZARDOUS SUBSTANCE. NON-DANGEROUS GOODS.
Hazards Overview Avoid contact with conductive gels during electronic measurements e.g. EEG or ECG. Do not use on patients undergoing MRI (Magnetic Resonance Imaging) examinations. Contains silver metal. Possible sensitisation in some individuals.
Risk Phrases Not applicable.
Safety Phrases Not applicable.

SECTION 3 – COMPOSITION / INFORMATION ON INGREDIENTS

Chemical Name	CAS No.
<i>Coating on dressing :</i>	
Silver	7440-22-4
Silver Oxide	20667-12-3
<i>Acticoat / Acticoat 7 fabric consists of:</i>	
High Density Polyethylene Mesh	N/A
Rayon / Polyester Gauze	N/A

MATERIAL SAFETY DATA SHEET
ACTICOAT PRODUCT RANGE
MSDS No. 179

<i>Acticoat Absorbent non-woven fabric containing:</i> Calcium Alginate	9005-35-00
<i>Acticoat Moisture Control / Acticoat Site fabric consists of:</i> Polyurethane wound contact layer	N/A
Central polyurethane Foam Layer	N/A
Blue pigmented polyurethane backing film	N/A

SECTION 4 – FIRST AID MEASURES

Swallowed	Obtain medical attention.
Eye	In the event of eye irritation, flush eyes with flowing water.
Skin	Obtain medical attention in the event that irritation develops or persists.
Inhaled	Not applicable.
First-Aid Facilities	Normal washroom facilities nearby.
Advice to Doctor	Treat symptomatically.

SECTION 5 – FIRE FIGHTING MEASURES

Hazards From Combustion Products	The product is a combustible solid which gives off toxic fumes, mainly carbon monoxide and carbon dioxide when ignited.
Suitable Extinguishing Media	According to the surrounding area. Any standard fire extinguisher may be used.
Fire-fighting Procedures	Wear self-contained breathing apparatus where large quantities of products are burning. Risk of explosion expected to be low.
Hazchem Code	Not applicable.

SECTION 6 – ACCIDENTAL RELEASE MEASURES

Spill / Leak Procedures	Treat used dressings as medical waste.
--------------------------------	--

MATERIAL SAFETY DATA SHEET**ACTICOAT PRODUCT RANGE****MSDS No. 179**

SECTION 7 – HANDLING AND STORAGE

Handling Advice Refer to instructions for use provided with product.
Storage Advice Store in cool (<25°C) dry environment, away from strong oxidising agents.
Keep away from moisture. Protect from light. Avoid extreme heat.

SECTION 8 – EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Standards ACGIH TLV (2000): Silver 0.1 mg/m³; Silver oxide 0.1 mg/m³.
LC₅₀ / LD₅₀: Silver – not available; Silver oxide Oral LD₅₀ Rat: 2820 mg/kg.
Engineering Controls Not required.
Personal Protection No special handling procedures required.

SECTION 9 – PHYSICAL AND CHEMICAL PROPERTIES

Appearance and Odour Solid fabric wound dressing.
Acticoat, Acticoat 7, Acticoat Moisture Control and Acticoat Site are blue on top, metallic silver on bottom.
Acticoat Absorbent is grey/brown in colour. No odour; metallic taste.

Melting Point / Boiling Point Melting / Sublimation Point: 120°C.
Vapour Pressure Not applicable.
Vapour Density (Air=1) Not applicable.
Specific Gravity >1
Flashpoint Not available.
Flammability Limits (FL) (%) Not available.
Autoignition Temperature Not available.
Solubility in Water (g/L) Silver coating is sparingly soluble in water.
For Acticoat Absorbent : fibre will gel or dissolve in an aqueous solution containing sodium or certain ions.

Percent volatile by volume Less than 1%.

MATERIAL SAFETY DATA SHEET**ACTICOAT PRODUCT RANGE****MSDS No. 179**

SECTION 10 – STABILITY AND REACTIVITY

Stability	Stable.
Conditions to Avoid	Naked flames and ignition sources.
Incompatible Materials	The silver coating is incompatible with these products: acetylene, acids, amines, ammonia, hydrazine, metal sulphides, peroxides, strong oxidisers.
Hazardous Decomposition Products	Carbon monoxide, carbon dioxide.
Hazardous Polymerisation	Not expected to occur.

SECTION 11 - TOXICOLOGICAL INFORMATION

Acute Effects	
Swallowed	Ingestion of the product may cause a choking hazard or internal blockage.
Eye	Corneal irritation may occur as a result of mechanical abrasion.
Skin	Skin irritation may result as an allergic reaction to the silver content of the dressing. Some individuals may be sensitised to silver metal, resulting in temporary discolouration of affected areas or rashes and irritation.
Inhaled	The product does not present an inhalation hazard.
Chronic Effects	Excessive or chronic exposure to silver by skin contact or ingestion may cause Argyrosis – a local or generalised discolouration of the eyes, skin and mucous membranes, where grey or blue patches of pigmentation are formed.

SECTION 12 – ECOLOGICAL INFORMATION

Environmental Issues	No environmental issues expected. Outer cardboard packaging material may be separated for recycling.
Persistence & Degradability	Not available.
Mobility	Not available.

SECTION 13 – DISPOSAL CONSIDERATIONS

Method of Disposal	In accordance with local waste regulations. Used material must be disposed of in accordance with infectious waste regulations.
---------------------------	--

MATERIAL SAFETY DATA SHEET**ACTICOAT PRODUCT RANGE****MSDS No. 179**

SECTION 14 – TRANSPORT INFORMATION

Classification	Not classified as a Dangerous Good according to the Australian Code for Transport of Dangerous Goods by Road and Rail.
UN No.	Not applicable.
Class	Not applicable.
Subsidiary Risk(s)	Not applicable.
Packing Group	Not applicable.
Hazchem Code	Not applicable.
Shipping Name	Not applicable.
Other	Not applicable.

SECTION 15 – REGULATORY INFORMATION

Poison Schedule Not a Scheduled Poison.**SECTION 16 – OTHER INFORMATION**

CONTACT POINT :	CUSTOMER SERVICE	Telephone: (03) 8540 6777
Date of Issue:	13/06/2008	
<p>This MSDS summarises our best knowledge of the health and safety hazard information of the product and how to safely handle and use the product in the workplace. Each user should read this MSDS and consider the information in the context of how the product will be handled and used in the workplace including in conjunction with other products. If clarification or further information is needed to ensure that an appropriate risk assessment can be made, the user should contact this company. Our responsibility for products sold is subject to our standard terms and conditions, a copy of which is sent to our customers and is also available on request.</p>		

END OF DOCUMENT

Item #	Cat #	Abbrev	MSDS Code	Description 1	Description 2
426218	20141	SNWOND		2 DRESSING, ACTICOAT 7DAY 4X5	(5/PK 6PK/CS) SNWOND
426222	20181	SNWOND		2 ACTICOAT, 3/4"X12" ROPE (5PK 1	2PK/CS) SNWOND
426232	20211	SNWOND		2 ACTICOAT, 5"X5" (5PK 4PK/CS)	SNWOND
426244	20381	SNWOND		2 ACTICOAT, 4"X5" (5PK 12PK/CS)	SNWOND
439889	20241	SNWOND		2 DRESSING, ACTICOAT 7DAY 6X6 (5	/PK 6PK/CS) SNWOND
440433	20101	SNWOND		2 DRESSING, ACTICOAT 4"X4" (12/P	K 4PK/CS) SNWOND
440434	20201	SNWOND		2 ACTICOAT, 4"X8" (12/PK 4PK/CS)	SNWOND
440435	20301	SNWOND		2 ACTICOAT, 8"X16" (6/PK 4PK/CS)	SNWOND
440436	20401	SNWOND		2 ACTICOAT, 16"X16" (6/PK 2PK/CS) SNWOND
440437	20501	SNWOND		2 ACTICOAT, 4"X48" (6/PK 2PK/CS)	SNWOND
473239	20211	SNWOND		2 DRESSING, ACTICOAT 5"X5" (5/BX) SNWOND
517026	20151	SNWOND		2 DRESSING, ACTICOAT W/NANOCRYST	ALLINE 5X5 (5BX/CS) SNWOND
517028	20341	SNWOND		2 DRESSING, ACTICOAT W/NANOCRYST	STALLINE 2X2 (20BX/C SNWOND
517031	20601	SNWOND		2 DRESSING, ACTICOAT W/NANOCRY	2"X2" (20BX/CS) SNWOND
535094	20111	SNWOND		2 DRESSING, ACTICOAT MOIST CNTRL	2X2 (10/BX) SNWOND
535095	20311	SNWOND		2 DRESSING, ACTICOAT MOIST CNTRL	4X4 (10/BX) SNWOND
585654	20611	SNWOND		2 DRESSING, ACTICOAT SITE DISC	SNWOND
632506	66021771	SNWOND		2 DRESSING, ACTICOAT POST-OP	SNWOND
632975	20611	SNWOND		2 DRESSING, ACTICOAT SITE RND 1"	(10/BX 10BX/CS) D/S SNWOND
634332	66021770	SNWOND		2 DRESSING, ACTICOAT PST-OP 4"X	5" (5/BX) SNWOND
634333	66021772	SNWOND		2 DRESSING, ACTICOAT PST-OP 4"X1	0" (5/BX) SNWOND
634397	66021771	SNWOND		2 DRESSING, ACTICOAT PST-OP 4"X	8" (5/BX) SNWOND
641854	66021771	SNWOND		2 DRESSING, ACTICOAT POST-OP (5/	BX 5BX/CS) SNWOND
641855	66021772	SNWOND		2 DRESSING, ACTICOAT PST-OP 4"X1	0" (5/BX 5BX/CS) SNWOND
641856	66021770	SNWOND		2 DRESSING, ACTICOAT PST-OP 4"X4	3/4" (5/BX 10BX/CS) SNWOND
697286	20381	SNWOND		2 DRESSING, ACTICOAT ABSRBNT 4X5	(5/BX 12BX/CS) SNWOND
702280	20241	SNWOND		2 ACTICOAT, 6"X6" (5/BX 4BX/CS)	(STERLING) SNWOND
704963	66800406	SNWOND		2 ACTICOAT, FLEX 3 4"X4" (12/BX)	D/S SNWOND
710957	66800405	SNWOND		2 DRESSING, ACTICOAT FLEX 7 4X5	(5/BX) SNWOND
711288	66800406	SNWOND		2 DRESSING, ACTICOAT FLEX 3 ANTI	MICRO 4X4 (12/BX 4BX SNWOND
711289	66800402	SNWOND		2 ACTICOAT, FLEX 3 ANTIMICRO DRS	2X2 (5/BX 20BX/CS) SNWOND
711290	66800417	SNWOND		2 ACTICOAT, FLEX 3 ANTIMICRO DRS	4X8 (12/BX 4BX/CS) SNWOND
712167	66800418	SNWOND		2 DRESSING, ACTICOAT FLEX3 8X16"	(6/BX 4BX/CS) SNWOND
712168	66800433	SNWOND		2 DRESSING, ACTICOAT FLEX3 16X1	6" (6/BX 2BX/CS) SNWOND
712169	66800434	SNWOND		2 DRESSING, ACTICOAT FLEX3 4X48"	(6/BX 2BX/CS) SNWOND
712170	66800418	SNWOND		2 DRESSING, ACTICOAT FLEX3 8X16"	(6/BX 4BX/CS) SNWOND
719349	66800405	SNWOND		2 DRESSING, ACTICOAT FLEX 7 4X5	(5/BX 6BX/CS) SNWOND
727296	66800403	SNWOND		2 DRESSING, ACTICOAT FLEX 7 ANTI	MICRO 2X2 (5/BX 20BX SNWOND
727742	66800544	SNWOND		2 DRESSING, ACTICOAT FLEX 7 RIBB	ON 1X24 (5/BX 10BX/C SNWOND
763809	66800427	SNWOND		2 DRESSING, ACTICOAT FLEX 7 ANTI	M 6X6 (5/BX 4BX/CS) SNWOND